

Sterile AROSuture™ by Needle Type, 2-0 to 7-0

AROSurgical
INSTRUMENTS CORPORATION


220 Newport Center Drive, Suite 11101

Newport Beach, CA 92660


Phone: 949-640-9932 or 888.430.7888

Fax: 949.640.9934 or 888.430.7889

www.arosurgical.com




	STERILE 2-0 - 5-0	SIZE	N.LENGTH	CIRCLE / CURVE	MATERIAL	S.LENGTH
IN STOCK	SP27A02N-45	2-0	27mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	SP19A03N-45	3-0	19mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	SP24A03N-45	3-0	24mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	SP15A04H-45	4-0	15mm	1/2circle 180°	Black Mono	45cm/18in
IN STOCK	SP15A04N-45	4-0	15mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	SP19A04N-45	4-0	19mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	SP13A05H-45	5-0	13mm	1/2circle 180°	Black Mono	45cm/18in
IN STOCK	SP13A05N-45	5-0	13mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	SP17A05N-45	5-0	17mm	3/8circle 135°	Black Mono	45cm/18in



	STERILE 4-0 to 7-0	SIZE	N.LENGTH	CIRCLE / CURVE	MATERIAL	S.LENGTH
IN STOCK	CC13A04N-45	4-0	13mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	CC13A05N-45	5-0	13mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	CC17A05N-45	5-0	17mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	CC11A06N-45	6-0	11mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	CC13A06N-45	6-0	13mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	CC11A07N-45	7-0	11mm	3/8circle 135°	Black Mono	45cm/18in



	STERILE 4-0 to 6-0	SIZE	N.LENGTH	CIRCLE / CURVE	MATERIAL	S.LENGTH
IN STOCK	E15A04L-45	4-0	15mm	7/16circle 157.5°	Black Mono	45cm/18in
IN STOCK	E15A04N-45	4-0	15mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	E19A04N-45	4-0	19mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	E13A05L-45	5-0	13mm	7/16circle 157.5°	Black Mono	45cm/18in
IN STOCK	E13A05N-45	5-0	13mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	E19A05N-45	5-0	19mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	E11A06N-45	6-0	11mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	E13A06N-25	6-0	13mm	3/8circle 135°	Black Mono	25cm/10in
IN STOCK	E17A06N-45	6-0	17mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	E19A06N-45	6-0	19mm	3/8circle 135°	Black Mono	45cm/18in



	STERILE 4-0 to 7-0	SIZE	N.LENGTH	CIRCLE / CURVE	MATERIAL	S.LENGTH
IN STOCK	BP13A04N-45	4-0	13mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	BP11A05N-45	5-0	11mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	BP13A05N-45	5-0	13mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	BP11A06N-45	6-0	11mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	BP13A06N-45	6-0	13mm	3/8circle 135°	Black Mono	45cm/18in
IN STOCK	BP13A07N-45	7-0	13mm	3/8circle 135°	Black Mono	45cm/18in